

SINGLE USE



D3010

Our improved mirror design features a zinc alloy tip for percussion testing, making the D3010 ideal for extraction or root analysis. In addition to our Front Face Mirror is our Sickle/CPITN-C Probe and our London College Tweezers.

Please see website link for kit details [here](#).

Single Use Dental Instruments Pty Ltd

1/24 Hoopers Road, Kunda Park
QLD 4556 Australia
sales@sudinstruments.com
07 5445 6486

ERGONOMIC • PRECISE • QUALITY • STERILE • SINGLE USE

FF53

Front Face Mirror


Description:

An essential item in any dental practice our Front Face Mirror is new and improved. With a focus on comfort, our ergonomic handle features additional ridges and grooves to improve grip stability. With the addition of a zinc alloy tip for percussion testing and an anti-fog mirror surface.

Features:

- Zinc alloy tip for percussion testing
- Ergonomic handle with ridges and grooves for improved grip
- Latex-free
- Arrives sterile and ready for use
- GS1 Data Matrix
- CE and TGA Certified

Materials: Zinc Alloy, Acrylonitrile butadiene styrene (ABS), Polycarbonate (PC), Aluminiumbutadiene styrene (ABS)



PB85

Sickle/CPITN-C Probe


Description:

Our PB85 combines an explorer and CPITN-C probe into one instrument. Featuring precision edges and easily identifiable measurements, this probe is a convenient addition to diagnostic assessments.

Features:

- Combination instrument
- Precision manufactured stainless steel tips
- Ergonomic handle with ridges and grooves for improved grip
- Calibrations; CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



TW82

London College Tweezers


Description:

Strong and sturdy with enhanced grip features on the handle and tips. Our premium London College Tweezers are an essential during procedures.

Features:

- Ergonomic handle
- Serrated non-slip tips
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel



Our Products

Our unique, single-use instruments were designed with simplicity in mind. Featuring a broad variety of diagnostic and restorative instruments, our goal is to remove daily stress for clinicians and improve patient safety. Not only convenient, we have created a product range that is focused on quality and sustainability. Our single use instruments arrive sterile and ready for immediate use. The ergonomic handle is made from ABS plastic and features ridges and grooves for an improved grip. The stainless-steel stem provides additional weight and durability to ensure it is fit for purpose. Precision manufactured stainless steel tips replicate the design and function of re-usable instruments.

Packaging and materials

We are committed to reducing the environmental impact of our products. As such our products and packaging are manufactured with recyclable materials. Additionally, we source recycled cardboard to create our benchtop dispensers and shipping containers.

Quality assurance

All our instrument manufacturing has CE, TGA and ISO certification to ensure quality compliance with industry standards.

Tyvek® medical dialysing strip

A Tyvek® medical dialysing strip is incorporated into all our product sleeves. Comprising of flashspun and bonded high-density polyethylene (HDPE) filaments, the strip provides a microbial barrier and exceptional tear and puncture resistance to maintain product integrity during transit.

