

SINGLE USE



S4001

One of our best-sellers, this kit features our unique combination Flat Plastic/Ball Burnisher, giving you two instruments in one in addition to our Small Spoon Excavator (1.2mm), Medium Spoon Excavator (1.8mm) and Ladmore Plugger. This kit is great value, the S4001 is a go-to for your restorative needs.

Single Use Dental Instruments Pty Ltd

1/24 Hoopers Road, Kunda Park
QLD 4556 Australia
sales@sudinstruments.com
07 5445 6486

Please see website link for kit details [here](#).

ERGONOMIC • PRECISE • QUALITY • STERILE • SINGLE USE

XC67

Small Spoon Excavator (1.2mm)


Description:

One of our bestsellers, our Small Spoon Excavator is a popular choice for removing caries within deep cavities.

Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Small spoon diameter; 1.2mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- CS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



XC58

Medium Spoon Excavator (1.8mm)


Description:

Gone are the days of dull and blunt instruments. Our single-use XC58 will have sharp edges each and every time. Precision made, our medium spoon excavator adds convenience to your procedures.

Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Medium spoon diameter; 1.8mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



CB59

Flat Plastic/ Ball Burnisher


Description:

The combined flat plastic and ball burnisher is brilliant for placing and blending composite material.

Features:

- Combination instrument
- Precision manufactured stainless steel tips
- Ergonomic handles with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



PL57

Ladmore Plugger


Description:

Our Ladmore Plugger features both medium and large tips which are ideal for placing and adapting composite.

Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Tip diameters; 1.0mm/2.1mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- CS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



Our Products

Our unique, single-use instruments were designed with simplicity in mind. Featuring a broad variety of diagnostic and restorative instruments, our goal is to remove daily stress for clinicians and improve patient safety. Not only convenient, we have created a product range that is focused on quality and sustainability. Our single use instruments arrive sterile and ready for immediate use. The ergonomic handle is made from ABS plastic and features ridges and grooves for an improved grip. The stainless-steel stem provides additional weight and durability to ensure it is fit for purpose. Precision manufactured stainless steel tips replicate the design and function of re-usable instruments.

Packaging and materials

We are committed to reducing the environmental impact of our products. As such our products and packaging are manufactured with recyclable materials. Additionally, we source recycled cardboard to create our benchtop dispensers and shipping containers.

Quality assurance

All our instrument manufacturing has CE, TGA and ISO certification to ensure quality compliance with industry standards.

Tyvek® medical dialysing strip

A Tyvek® medical dialysing strip is incorporated into all our product sleeves. Comprising of flashspun and bonded high-density polyethylene (HDPE) filaments, the strip provides a microbial barrier and exceptional tear and puncture resistance to maintain product integrity during transit.

