

# SINGLE USE



## S6001

For complex restorations, our S6001 has everything you need. Reducing instrument preparation into one easy step, this kit allows you to spend your time focused directly on the patient. Inclusive of the Six Piece Restorative Kit is our Ball Burnisher, Carver, Medium Spoon Excavator (1.8mm), Flat Plastic, Amalgam Plugger (Large/Small) and Ladmore Plugger which features both medium and large tips.

Please see website link for kit details [here](#).

Single Use Dental Instruments Pty Ltd

1/24 Hoopers Road, Kunda Park  
QLD 4556 Australia  
sales@sudinstruments.com  
07 5445 6486

ERGONOMIC • PRECISE • QUALITY • STERILE • SINGLE USE

# BB56

## Ball Burnisher


### Description:

The perfect blending tool. Our ball burnisher has both small and large tip ends to effectively condense and shape the composite against the cavity wall surfaces.

### Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- CS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



# CV55

## Carver


### Description:

A brand new instrument for every patient, our CV55 is sharp and precise.

### Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



# XC58

## Medium Spoon Excavator (1.8mm)


### Description:

Gone are the days of dull and blunt instruments. Our single-use XC58 will have sharp edges each and every time. Precision made, our medium spoon excavator adds convenience to your procedures.

### Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Medium spoon diameter; 1.8mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



# FP54

## Flat Plastic


### Description:

Our double ended Flat Plastic features long, flat tips that are used to carry and shape composite materials.

### Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Tip diameters; 2.2mm/1.8mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- CS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



# PL61

## Amalgam Plugger (Small/Large)


### Description:

The single-use Amalgam Plugger has both a large and a small plugging surface. The combination is useful for both broad and deep amalgam condensing to ensure a high-quality filling.

### Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Tip diameters; 2.5mm/1.5mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



# PL57

## Ladmore Plugger


### Description:

Our Ladmore Plugger features both medium and large tips which are ideal for placing and adapting composite.

### Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Tip diameters; 1.0mm/2.1mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- CS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



# Our Products

Our unique, single-use instruments were designed with simplicity in mind. Featuring a broad variety of diagnostic and restorative instruments, our goal is to remove daily stress for clinicians and improve patient safety. Not only convenient, we have created a product range that is focused on quality and sustainability. Our single use instruments arrive sterile and ready for immediate use. The ergonomic handle is made from ABS plastic and features ridges and grooves for an improved grip. The stainless-steel stem provides additional weight and durability to ensure it is fit for purpose. Precision manufactured stainless steel tips replicate the design and function of re-usable instruments.

## Packaging and materials

We are committed to reducing the environmental impact of our products. As such our products and packaging are manufactured with recyclable materials. Additionally, we source recycled cardboard to create our benchtop dispensers and shipping containers.

## Quality assurance

All our instrument manufacturing has CE, TGA and ISO certification to ensure quality compliance with industry standards.

## Tyvek® medical dialysing strip

A Tyvek® medical dialysing strip is incorporated into all our product sleeves. Comprising of flashspun and bonded high-density polyethylene (HDPE) filaments, the strip provides a microbial barrier and exceptional tear and puncture resistance to maintain product integrity during transit.

