



SINGLE USE

DENTAL INSTRUMENTS

By listening to customers and continuously working to improve our offering, our range has evolved to feature a wide variety of instruments and kits ideal for every procedure.

Individual Instruments

Using individually wrapped instruments is a great way to transition your clinic towards implementing single-use solutions. Only open those you need for procedures and take minimising waste to a whole new level.

Restorative Kits

Our restorative instrument kits will leave you staring at the autoclave wondering why you spend so much money maintaining it every year. Using our thoughtfully designed combinations alleviates the pressure and expense of reprocessing. Remember if you can't find a restorative combination that works for you, we'll make you one. That simple!

Diagnostic Kits

Have you been entertaining the idea of single-use instruments? Equipping your clinic with our readily available diagnostic kits allows your nurse to spend less time setting up and more time chair side, focusing on the patient. Who would complain about that?

Individual Instrument Range	4–19
Burnishers	6
Carvers	7
Excavators	8
Gracey Curettes	9
Universal Curettes	10
Flat Plastics	11
Mirror	12
Pluggers	13
Probes	14–16
Scalers	17
Tweezers	18
Other	19
Diagnostic Instrument Range	20–31
Two Piece Diagnostic Kits	22–27
Three Piece Diagnostic Kits	28–32
Complete Periodontal Diagnostic Kit	33
Restorative Instrument Range	35–49
Two Peice Restorative Kits	35–39
Three Peice Restorative Kits	40–41
Four Peice Restorative Kits	42–45
Five Peice Restorative Kits	47
Six Peice Restorative Kits	48–49

INDIVIDUAL

INSTRUMENT RANGE

Burnishers

Carvers

Excavators

Curettes

Flat Plastics

Mirrors


Pluggers




Probes

Scalers

Tweezers

Other





BB49

BEST SELLER

PFI49 Burnisher

The PFI49 burnisher is ideal for placing and contouring restorations. Having both a rounded and pointed tip makes this instrument useful for placing composite in small increments prior to curing.

BB49 | 20 Per Box

BB56

BEST SELLER

Ball Burnisher

The perfect blending tool. Our ball burnisher has both small and large tip ends to effectively condense and shape the composite against the cavity wall surfaces.

BB56 | 20 Per Box

CB59

Flat Plastic & Ball Burnisher

The combined flat plastic and ball burnisher is brilliant for placing and blending composite material.

CB59 | 20 Per Box

CV55

Carver

Being single-use, our carver tips are always sharp and perfectly pointed.

CV55 | 20 Per Box

CV66

½ Hollenback

The CV66 is a versatile instrument which is very useful for placing, carving and contouring amalgam.

CV66 | 20 Per Box



CX66

BEST SELLER

**Extra Small
Spoon Excavator**

The go-to instrument for clearing out the smallest micro-dental reparations.

CX66 | 20 Per Box

XC67

BEST SELLER

**Small
Spoon Excavator**

One of our best sellers, our small spoon excavator is the top choice when you want to have extra control over caries removal in deep cavities.

XC67 | 20 Per Box

XC58

BEST SELLER

**Medium
Spoon Excavator**

Gone are the days of dull and blunt instruments. Precision made, our single-use XC58 will have sharp edges each and every time.

XC58 | 20 Per Box

XC68

**Large
Spoon Excavator**

The large curved surface on our XC68 excavator is ideal for removing broad-scale decay.

XC68 | 20 Per Box

GC11

Gracey Curette | 11 & 12

Using disposable curettes gives you the optimum working end in each appointment. Our GC11 is a great choice for finishing subgingival calculus removal.

GC11 | 20 Per Box

GC13

Gracey Curette | 13 & 14

Adopting the use of our single-use GC13 means the lower cutting edge will be prime for scaling, every time. The smooth edges on our curettes are perfect for removing subgingival calculus and getting into pockets without harming the patient.


GC13 | 20 Per Box

GC17

Gracey Curette | 17 & 18

Our disposable GC17 will always have the perfect angle for reaching distal surfaces of posterior teeth. The thick handle design provides you with comfort and accentuated grip throughout procedures.

GC17 | 20 Per Box

**UC02****BEST SELLER****Universal Curette | 2R & 2L**

Never have a dull cutting edge again when you make the switch to our single-use UC02. Fresh out the pack, they provide you with reliable quality in your scaling procedures.

UC02 | 20 Per Box

UC04**Universal Curette | 4R & 4L**

Sterile and ready for use, our UC04 is the convenient choice for incisors, canines and premolars.

UC04 | 20 Per Box

FP48**Small Flat Plastic**

About half the size of our original flat plastic, the tips on our FP48 are ideal for small and narrow cavities.

FP48 | 20 Per Box

FP54**BEST SELLER****Flat Plastic**

Our double-ended flat plastic features long, flat tips that are helpful for carrying and shaping composite materials.

FP54 | 20 Per Box

CB59**Flat Plastic & Ball Burnisher**

The combined flat plastic and ball burnisher is brilliant for placing and blending composite material.

CB59 | 20 Per Box




FF53

BEST SELLER

Front Face Mirror with Percussion Tip

An essential item in any dental practice, our front face mirror is new and improved. With a focus on comfort, our ergonomic handle features additional ridges and grooves to improve grip stability. With the addition of a zinc alloy tip for percussion testing and an anti-fog mirror surface... What more could you ask for?

FF53 | 20 Instruments Per Box



PL57

Ladmore Plugger

Our ladmore plugger features both medium and large tips which are ideal for placing and adapting composite.

PL57 | 20 Per Box

PL61

Amalgam Plugger

The single-use amalgam plugger has a large and a small plugging surface. The combination is useful for both broad and deep amalgam condensing to ensure a strong filling.

PL61 | 20 Per Box



PB84

**Explorer No. 8
& CPITN-C Probe**

Ensuring that our product range caters for everyone our PB84 is popular amongst dentists trained within the UK. CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

PB84 | 20 Per Box

PB85

**Sickle Explorer
& CPITN-C Probe**

Our PB85 combines an explorer and CPITN-C probe into one instrument. Featuring precision edges and easily identifiable measurements, this probe is a convenient addition to diagnostic assessments.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

PB85 | 20 Per Box

PB89

BEST SELLER

**Sickle Explorer
& Perio Probe**

A customer favourite, our perio probe & sickle has easily identifiable markings to assist in delivering accurate assessments. This combination instrument can be purchased individually or within our top selling three-piece diagnostic kits: D3001 and D3006.

Perio Probe: 3/5/7/9/11/13mm.

PB89 | 20 Per Box

PB90

**Williams
& CPITN-C Probe**

The PB90 features a Williams tip with the 4mm and 6mm markings omitted for easy identification and a CPITN-C probe with additional 8.5–11.5mm markings.

Williams: 1/2/3/5/7/8/9/10mm.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

PB90 | 20 Per Box



PB91

UNC-15 & CPITN-C Probe

Combining the popular ball tip/WHO probe with our UNC-15, this probe is an excellent choice when measuring pocket depth.

UNC-15: 1/2/3/4/5/6/7/8/9/10/11/12/13/14/15.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

PB91 | 20 Per Box




PB93

Explorer No. 6

Our single-ended No. 6 Probe is useful for many purposes. Whether you're inspecting filling margins or checking cavities prior to filling, this probe is in prime condition and ready to go.

PB93 | 20 Per Box



SC02

Sickle Scaler 204S

A posterior sickle scaler designed to remove supragingival interproximal calculus. Newly improved, our SC02 features a thick handle with ridges and grooves for improved comfort and grip.

SC02 | 20 Per Box



SC03

BEST SELLER

Hygienist Scaler | H6 & H7

A favourite amongst hygienists this scaler is designed for interproximal cleaning between anterior teeth and premolars. Featuring a stainless-steel stem and ergonomic handle to ensure stability and control during procedures.

SC03 | 20 Per Box



TW81

Examination Tweezers

Our economy tweezer is a popular choice for both diagnostic and restorative consultations. Featuring serrated, non-slip tips to ensure you have complete control during procedures.

TW81 | 20 Per Box




TW82

BEST SELLER

London College Tweezers

Strong and sturdy with enhanced grip features on the handle and tips. Our premium London College Tweezers are an essential tool in most procedures.

TW82 | 20 Per Box



DA69

BEST SELLER

Dycal Applicator

Our single-ended dycal applicator is strong and slender to ensure accurate placement of lining materials.

DA69 | 20 Instruments Per Box



SP65

BEST SELLER

Spatula

The single-use spatula is slim, sleek and refined. Perfectly clean every time, these spatulas remove the reprocessing hassle after mixing and spreading sticky dental materials.

SP65 | 20 Instruments Per Box

DIAGNOSTIC

INSTRUMENT RANGE

Two Piece Kits

Three Piece Kits

Complete Diagnostic Kit





D2004

Two Piece | Diagnostic Kit

Great for complete assessments. Our combination probe has a CPITN tip ideal for basic periodontal exams and a Williams for in-depth measurements. The probe markings are easily identifiable and in accordance with WHO guidelines.

- **Front Face Mirror with Percussion Tip**
- **Williams & CPITN-C Probe**

D2004 | 150 Kits Per Box



D2006

BEST SELLER

Two Piece | Diagnostic Kit

One of our best sellers, the combination perio probe & sickle is beneficial for both periodontal screening and caries detection.

- **Front Face Mirror with Percussion Tip**
- **Sickle Explorer & Perio Probe**

D2006 | 150 Kits Per Box





D2009

BEST SELLER

Two Piece | Diagnostic Kit

Everything you need for a standard examination. New and improved, our front face mirror allows for high-quality illumination and retraction. Precision manufactured with stainless steel tips, our CPITN-C probe features easily identifiable markings in accordance with WHO guidelines.

- **Front Face Mirror with Percussion Tip**
- **Sickle & CPITN-C Probe**

D2009 | 150 Kits Per Box



D2014

BEST SELLER

Two Piece | Diagnostic Kit

One of our best sellers, the combination perio probe & sickle is beneficial for both periodontal screening and caries detection.

- **Front Face Mirror with Percussion Tip**
- **Explorer No. 8 & CPITN-C Probe**

D2014 | 150 Kits Per Box





D2015

Two Piece | Diagnostic Kit

Combining a Front Face Mirror with CPITN-C probe & UNC-15, the D2015 kit is designed for clear visibility and accurate periodontal screening.

- **Front Face Mirror with Percussion Tip**
- **UNC-15 & CPITN-C Probe**

D2015 | 150 Kits Per Box



D2016

Two Piece | Diagnostic Kit

Featuring a Front Face Mirror and UNC-15 Sickie Explorer, the D2016 kit is ideal for precise detection of caries and calculus during examinations.

- **Front Face Mirror with Percussion Tip**
- **Sickle Explorer & UNC-15**

D2016 | 150 Kits Per Box





D3001

BEST SELLER

Three Piece | Premium Diagnostic Kit

Our D3001 is flying off the shelves for good reason. Featuring our London college tweezers, this three-piece combination assists in delivering premium quality examinations.

- **Front Face Mirror with Percussion Tip**
- **Sickle Explorer & Perio Probe**
- **London College Tweezers**

D3001 | 100 Kits Per Box




D3003

Three Piece | Orthodontic Kit

Our orthodontic exam kits will be your go-to for quick assessments. The small ruler size is ideal for both paediatric and adult consultations.

- **Front Face Mirror with Percussion Tip**
- **Sickle Explorer & Perio Probe**
- **Orthodontic Ruler**

D3003 | 100 Kits Per Box





D3006

Three Piece | Diagnostic Kit

With a new and improved mirror and examination tweezer you can be confident in the quality of our D3006.

- Front Face Mirror with Percussion Tip
- Sickle Explorer & Perio Probe
- Examination Tweezers

D3006 | 100 Kits Per Box



D3010

Three Piece | Premium Diagnostic Kit

Our improved mirror design features a zinc alloy tip for percussion testing, making the D3010 ideal for indicating extraction or root treatment.

- Front Face Mirror with Percussion Tip
- Sickle Explorer & CPITN-C Probe
- London College Tweezers

D3010 | 100 Kits Per Box





D3014

Three Piece | Diagnostic Kit

Featuring our new and improved front face mirror and examination tweezers as well as a CPITN-C ball tip probe.

- Front Face Mirror with Percussion Tip
- Explorer No. 8 & CPITN-C Probe
- Examination Tweezers

D3014 | 100 Kits Per Box



D3550

BEST SELLER

Complete Periodontal | Diagnostic Kit

Convenient for mobile dentists and domiciliary care. With everything at your fingertips, our diagnostic pack is making patient examinations quick and easy.

- Front Face Mirror with Percussion Tip
- Sickle Explorer & Perio Probe
- Tweezers
- 2 x Face Mask with Ear Loops
- Patient Bib with Ties
- 2 x Cotton Roll No.2
- Tissue wipe
- Tray | 19.5cm x 12.5cm x 2cm

D3550 | 25 Kits Per Box

RESTORATIVE

INSTRUMENT RANGE

Two Piece Kits

Three Piece Kits

Four Piece Kits

Five Piece Kits

Six Piece Kits





S2000

Two Piece | Restorative Kit

Suitable for any filling, this two-piece kit is a popular choice.

- **Ball Burnisher**
- **Flat Plastic**

S2000 | 150 Kits Per Box



S2002

BEST SELLER

Two Piece | Restorative Kit

Say goodbye to the days of scrubbing mixing materials from your instruments with our S2002. This single-use mixing kit gives you spotless tools every time.

- **Dycal Applicator**
- **Spatula**

S2002 | 150 Kits Per Box





S2003

Two Piece | Restorative Kit

This two-piece kit is popular for prepping and polishing a cavity restoration.

- **Ball Burnisher**
- **Small Spoon Excavator**

S2003 | 150 Kits Per Box



S2008

Two Piece | Restorative Kit

An ideal kit for placing and finishing Amalgam. Having a new plugger in every procedure means you never have to check for leftover residue after reprocessing. The small and large tips allow you to pack and contour different sized cavities.

- **Amalgam Plugger**
- **Carver**

S2008 | 150 Kits Per Box





S3000

BEST SELLER

Three Piece | Restorative Kit

One of our best sellers. With instruments for both placement and caries removal, this three-piece combination is perfect for most fillings.

- **Flat Plastic**
- **Ball Burnisher**
- **Medium Spoon Excavator**

S3000 | 100 Kits Per Box




S3012

Three Piece | Restorative Kit

From packing to perfection this kit has you covered. Your Amalgam restorations have never been so convenient.

- **Round Plugger**
- **Amalgam Plugger**
- **½ Hollenback**

S3012 | 100 Kits Per Box





S4000

Four Piece | Restorative Kit

This four-piece combination has everything you need to keep things clean and dry before placing composite. The inclusion of our premium London college tweezers makes this kit a favourite for restorative procedures.

- Flat Plastic
- Ball Burnisher
- Medium Spoon Excavator
- London College Tweezers

S4000 | 80 Kits Per Box



S4001

BEST SELLER

Four Piece | Restorative Kit

One of our best sellers this kit features our unique combination flat plastic & ball burnisher, giving you two instruments in one. A great value kit which is endorsed by dentists throughout Australia, it is a go-to for your restorative needs.

- Flat Plastic & Ball Burnisher
- Ladmore Plugger
- Small Spoon Excavator
- Medium Spoon Excavator

S4001 | 80 Kits Per Box





S4003

Four Piece | Restorative Kit

The S4003 just made prepping, placing and polishing a whole lot quicker.

- Flat Plastic
- Ball Burnisher
- Small Spoon Excavator
- Medium Spoon Excavator

S4003 | 80 Kits Per Box



S4004

Four Piece | Restorative Kit

With the inclusion of our small spoon excavator, the S4004 is perfect for paediatric restorations.

- Flat Plastic
- Ball Burnisher
- Small Spoon Excavator
- London College Tweezers

S4004 | 80 Kits Per Box





SHOP ONLINE NOW

singleuse.au



S5000

Five Piece | Restorative Kit

For a more complete Amalgam kit, look no further than our S5000.

- Flat Plastic
- Ball Burnisher
- Carver
- Amalgam Plugger | Small & Large
- Medium Spoon Excavator

S5000 | 60 Kits Per Box





S6000

Six Piece | Restorative Kit

The more complex the procedure the more instruments you require. Simplify your preparation process into one simple step. Our S6000 is ready when you are.

- Flat Plastic
- Ball Burnisher
- Burnisher
- 1/2 Hollenback
- Small Spoon Excavator
- Medium Spoon Excavator

S6000 | 50 Kits Per Box



S6001

Six Piece | Restorative Kit

For complex restorations, our S6001 has everything you need. Reducing instrument preparation into one easy step, this kit allows you to spend your time focused on the patient.

- Ladmore Plugger
- Amalgam Plugger | Small & Large
- Medium Spoon Excavator
- Flat Plastic
- Ball Burnisher
- Carver

S6001 | 50 Kits Per Box



SINGLE USE

07 5445 6486
sales@sudinstruments.com
singleusedentalinstruments.com

1/24 Hoopers Road
Kunda Park, QLD Australia 4556