


SINGLE USE



FF53

**Front Face Mirror
with Percussion Tip**

FF53 | Front Face Mirror with Percussion Tip



Description


An essential item in any dental practice our Front Face Mirror is new and improved. With a focus on comfort, our ergonomic handle features additional ridges and grooves to improve grip stability. With the addition of a zinc alloy tip for percussion testing and an anti-fog mirror surface.

Features

- Zinc alloy tip for percussion testing
- Ergonomic handle with ridges and grooves for improved grip
- Latex-free
- Arrives sterile and ready for use
- GS1 Data Matrix
- CE Certified

Materials

- Zinc Alloy
- Acrylonitrile butadiene styrene (ABS)
- Polycarbonate (PC)
- Aluminiumbutadiene styrene (ABS)



Our Products

Our unique, single-use instruments were designed with simplicity in mind. Featuring a broad variety of diagnostic and restorative instruments, our goal is to remove daily stress for clinicians and improve patient safety. Not only convenient, we have created a product range that is focused on quality and sustainability. Our single use instruments arrive sterile and ready for immediate use. The ergonomic handle is made from ABS plastic and features ridges and grooves for an improved grip. The stainless-steel stem provides additional weight and durability to ensure it is fit for purpose. Precision manufactured stainless steel tips replicate the design and function of re-usable instruments.

Packaging and materials

We are committed to reducing the environmental impact of our products. As such our products and packaging are manufactured with recyclable materials. Additionally, we source recycled cardboard to create our benchtop dispensers and shipping containers.

Quality assurance

All our instrument manufacturing has CE and ISO certification to ensure quality compliance with standards.

Tyvek® medical dialysing strip

A Tyvek® medical dialysing strip is incorporated into all our product sleeves. Comprising of flashspun and bonded high-density polyethylene (HDPE) filaments, the strip provides a microbial barrier and exceptional tear and puncture resistance to maintain product integrity during transit.

SINGLE USE

Single Use Dental Instruments Pty Ltd

07 5445 6486 | singleuse.au
sales@sudinstruments.com
1/24 Hoopers Road
Kunda Park, QLD Australia 4556
